

Instruments	Capabilities and preferred use								Intended Use		
	Proteomics				Metabolomics				Discovery	Translation/Routine Service	Clinical
	Multi Cell		Single Cell		Multi Cell		Single Cell		Highly flexible focusing to maximize data output and target discovery. Typical sample #: tens	Method optimization focusing to validate analytical method and molecular targets. Typical sample #: hundreds	Defined Standard Operational Procedures focusing on robustness and reproducibility. Typical sample #: thousands
	Profiling Shotgun	Targeted	FACS based	MALDI Imaging	Profiling Shotgun	Targeted	FACS based	MALDI Imaging			
Bruker's timsTOF fleX with ESI and MALDI	x		x	x	x		x	x	x		
Fisher Orbitrap Fusion Tribrid	x		x		x		x		x	x	x
Thermo Fisher Q Exactive Plus	x		x		x		x		x	x	
Thermo Fisher Q Exactive Plus	x		x		x		x			x	x
Agilent's 6495 Triple Quadrupole LC/MS		x				x			x	x	
Agilent's 6475 Triple Quadrupole LC/MS		x				x				x	x
Agilent 7250 GCMS QTOF with EI and CI sources					x				x	x	